



## Copper trace cuvette test 0.01-1.0 mg/L Cu, 20 tests

**Product #:** LCK529  
**AED Price:** Contact Hach

### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

**Convenient handling with 50 mm rectangular cuvettes + Dosicaps**

**Highest accuracy for analysis of ultra low range concentrations**

---

### Specifications

Description:	Copper trace, 50 mm cuvette test
Digestion required:	Yes, for Total Copper
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.01 - 1.0 mg/L Cu
Method:	Bathocuproine Disulphonic Acid
Number of tests:	20
Parameter:	Copper, trace
Platform:	LCK
Special remarks:	For total Copper digestion use Crack Set LCW902
Storage conditions:	15 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)