



Chromium (III and VI) trace cuvette test 0.005-0.25 mg/L Cr

Product #:

LCS313

AED Price:

Contact Hach

Hazardous



Hach cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

Specifications

According to standard:	EN ISO 11083, DIN 38405-D24
Description:	Chromium VI and total
Digestion required:	yes, for total Chromium
Measuring range:	0.005 - 0.25 mg/L Cr (VI)
Measuring range (2):	0.005 - 0.25 mg/L Cr total
Method:	Diphenylcarbazide
Number of tests:	25
Parameter:	Chromium, trace
Platform:	LCK
Special remarks:	for the evaluation 50mm semi-micro cuvettes needed
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)