



Nickel trace cuvette test 0.05-1.0 mg/L Ni, 20 tests

Product #:

LCK537

Hazardous

AED Price:

Contact Hach



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Convenient handling with 50 mm rectangular cuvettes

Highest accuracy for analysis of ultra low range concentrations

Specifications

Description:	Nickel trace, 50 mm rectangle cuvette test
Digestion required:	Yes, for total Nickel
EPA compliant:	N/A
Instrument:	DR3900, DR6000
Measuring range:	0.05 - 1.0 mg/L Ni
Method:	Dimethylglyoxime
Number of tests:	20
Parameter:	Nickel, trace
Platform:	LCK
Special remarks:	For total Nickel digestion use Crack Set LCW902
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)