



Manganese Reagent Set 0.2 - 5 mg/L / 0.02 - 1.0 mg/L Mn

Product #:

LCW032

AED Price:

Contact Hach

Hazardous



Manganese reagent set 0.2 - 5 mg/L Mn (Round cuvette or 10 mm rectangle cuvette) / 0.02 - 1.0 mg/L Mn (50 mm rectangle cuvette)

For determination of Manganese by the Formaldoxime method.

Specifications

Measuring range:	0.2 - 5 mg/L Mn
Measuring range (2):	0.02 - 1 mg/L Mn (50 mm rectangular cuvette)
Method:	Formaldoxime
Number of tests:	50
Parameter:	Manganese
Storage conditions:	2 °C - 8 °C (keep refrigerated)