



## Chromic acid baths cuvette test 0.5-450 g/L $\text{CrO}_3$ , 25 tests

Product #:

LCK213

Hazardous

AED Price:

Contact Hach



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

---

### Specifications

Description:	Acidic Chromium plating baths, galvanic chromium baths
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.5 - 5.0 g/L $\text{CrO}_3$
Measuring range (2):	50 - 450 g/L $\text{CrO}_3$
Method:	Intrinsic Baths Colour
Number of tests:	25
Parameter:	Chromic acid
Platform:	LCK
Storage conditions:	15 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)