



Bitter units (MEBAK) cuvette test ≥ 2 BU, 25 tests

Product #:

LCK241

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

For use with your Hach DR6000 spectrophotometer. LCK241 offers you the ability to measure international Bitter Units with less hassle and no chemical waste. Additionally, these cuvettes provide high quality results free from potential external influencers. They are designed to simplify your brewing analysis, with accurate results, every time.

Features include:

Ready-to-use, high-precision pre-dosed cuvette tests

Reduction of errors due to bar-coded cuvettes

RFID technology to document the shelf life of the cuvettes and ensure traceability of your analyses

No reagent blank needed

Innovative cuvette and reagent technique

Specifications

According to standard:	MEBAK II
Instrument:	DR5000, DR6000
Measuring range:	≥ 2 Bitter units
Method:	Analogous MEBAK-Method
Number of tests:	25
Parameter:	Bitter units
Platform:	LCK
Standard method:	MEBAK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)