



## Menthol in distillate cuvette test 0.5-15 mg Menthol/100 mL

Product #:

LYW185

AED Price:

Contact Hach

Hazardous



---

### Specifications

|                     |                              |
|---------------------|------------------------------|
| Measuring range:    | 0.5 - 15 mg/100 mL Menthol   |
| Method:             | p-Dimethylaminobenzaldehyde  |
| Number of tests:    | 25                           |
| Parameter:          | Menthol                      |
| Platform:           | LCK                          |
| Storage conditions: | 2 - 8 °C (keep refrigerated) |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)