



## Zirconium cuvette test, 6-60 mg/L Zr, 25 tests

Product #:

LCK364

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

The first specific, easy to use rapid test which helps to evaluate the Zirconium concentration during the passivation process.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

**Helps adjusting the right amount of Zirconium within the process**

**Helps minimising chemistry consumption and costs in Zirconium based coating processes**

**Helps ensuring a consistent high product quality**

---

### Specifications

EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	6 - 60 mg/L
Method:	SurTec/Hach Method
Number of tests:	25
Parameter:	Zirconium
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)