



Water Hardness cuvette test 1-20 °dH, 25 tests

Product #: LCK327
AED Price: Contact Hach

LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

Description:	Water Hardness, Calcium / Magnesium
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	1 - 20 °dH
Measuring range (2):	5 - 100 mg/L Ca
Measuring range (3):	3 - 50 mg/L Mg
Method:	Metalphthalein
Number of tests:	25
Parameter:	Water hardness
Platform:	LCK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)