



## Residual Hardness cuvette test 0.02-0.6 °dH, 24 tests

Product #:

LCK427

Hazardous

AED Price:

Contact Hach



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

---

## Specifications

Description:	Water Hardness, Calcium/Magnesium
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.02 - 0.6 °dH
Measuring range (2):	0.1 - 2 mg/L Ca
Measuring range (3):	0.15 - 2 mg/L Mg
Method:	Metaphthalein
Number of tests:	24 (test includes zero solution)
Parameter:	Water hardness, residual
Platform:	LCK
Storage conditions:	15 - 25 °C

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)