



Silica Reagent Set, 0.01-0.8 mg/L SiO₂

Product #:

LCW028

Hazardous

AED Price:

Contact Hach



Or 0.005-0.4 mg/L Si.

For determination of Silica by the Molybdenum Blue method.

Specifications

Description:	Reagent set for Silica acid/Silicium
Measuring range:	0.01 - 0.8 mg/L SiO ₂
Measuring range (2):	0.005 - 0.4 mg/L Si
Method:	Molybdenum Blue
Number of tests:	50
Parameter:	Silica
Sample size:	50 mL
Special Features:	50 tests
Special remarks:	Needed: round cuvettes LCZ906 or 50 mm rectangular cuvette
Storage conditions:	15 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)