



Hydrazine reagent set 0.01 - 2.0 mg/L N_2H_4

Product #:

LCW025

Hazardous

AED Price:

Contact Hach



For determination of Hydrazine by the 4-Dimethylaminobenzaldehyde method.

Specifications

Description:	reagent set for Hydrazine
Measuring range:	0.01 - 2.0 mg/L N_2H_4
Method:	4-Dimethylaminobenzaldehyde
Number of tests:	60
Parameter:	Hydrazine
Storage conditions:	15 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)