



## COD cuvette test - ISO 15705, 0-150 mg/L O<sub>2</sub>

Product #:

LCI500

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

**Exact concentrations of the ISO approved formula**

**Extremely accurate**

**Pre-dosed & ready to use**

**No contact with toxic reagent**

**Disposal/recycling service**

---

### Specifications

According to standard:	ISO 15705
Description:	Chemical Oxygen Demand
Digestion required:	Yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0 - 150 mg/L O <sub>2</sub>
Method:	Dichromate
Number of tests:	24 tests, one blank
Parameter:	COD
Platform:	LCK
Standard method:	ISO 15705
Storage conditions:	15 - 25 °C (protect from light)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)