



## COD cuvette test 100-2000 mg/L O<sub>2</sub>, 25 tests

**Product #:**

**LCK514**

**AED Price:**

Contact Hach

**Hazardous**



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	ISO 6060-1989, DIN 38409-H41-H44
Description:	Chemical Oxygen Demand
Digestion required:	Yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	100 - 2000 mg/L O <sub>2</sub>
Method:	Dichromate
Number of tests:	25
Parameter:	COD
Platform:	LCK
Storage conditions:	15 - 25 °C (protect from light)
Truecal:	Yes

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)