



## COD cuvette test 100-600 mg/L O<sub>2</sub>, 25 tests

Product #:

LCK714

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, for maximum safety, easy handling and photometric evaluation, environmentally friendly, approved for consent limits

**Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer and Truecal calibration**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### Specifications

According to standard:	ISO 6060-1989, DIN 38409-H41-H44
Description:	Chemical Oxygen Demand
Digestion required:	yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	100 - 600 mg/L O <sub>2</sub>
Method:	Dichromate
Number of tests:	25
Parameter:	COD
Platform:	LCK
Storage conditions:	15 - 25 °C (protect from light)
Truecal:	Yes

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)