



Formaldehyde trace cuvette test 0.01-1.0 mg/L H₂CO

Product #: LCS325
AED Price: Contact Hach

Specifications

Digestion required:	yes
Instrument:	DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.01 - 1.0 mg/L H ₂ CO
Method:	Acetylacetone
Number of tests:	24 (test includes a zero solution)
Parameter:	Formaldehyde, trace
Platform:	LCK
Special remarks:	For the evaluation 50 mm semi-micro rectangle-cuvettes are needed.
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)