



## Cationic surfactants cuvette test 0.2-2.0 mg/L, 25 tests

Product #:

LCK331

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

#### Maximum safety for users, thanks to the closed system and low amounts of reagents

#### Convenient and error-free dosing of the reagents

#### Barcode label for automatic recognition in the photometer

#### Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

Description:	Cetyltrimethylammoniumbromide (CTAB)
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.2 - 2.0 mg/L
Method:	Bromophenol Blue
Number of tests:	25
Parameter:	Surfactants, cationic
Platform:	LCK
Storage conditions:	15 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)