



## Phenols cuvette test 0.05-5.0 mg/L, 24 tests

Product #:

LCK345

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

**Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### Specifications

EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.05 - 5 mg/L Phenols
Method:	4-Nitroaniline
Number of tests:	24 (test includes zero solution)
Parameter:	Phenols
Platform:	LCK
Storage conditions:	15 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)