



Phenols cuvette test 5 - 150 mg/L, 24 tests

Product #:

LCK346

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 6439-1990, DIN 38409 H16
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	5 - 150 mg/L Phenols
Method:	4-Aminoantipyrine
Number of tests:	24 (test includes zero solution)
Parameter:	Phenols
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)