



## TOC cuvette test (purging method) 3-30 mg/L C, 25 tests

Product #:

LCK385

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	EN 1484, DIN 38409-H3
Description:	Total Organic Carbon
Digestion required:	Yes
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	3 - 30 mg/L C
Method:	Purging Method, Persulphate Digestion
Number of tests:	25
Parameter:	TOC
Platform:	LCK
Standard method:	EN 1484
Storage conditions:	2 - 8 °C (keep refrigerated)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)