



## AOX cuvette test 0.05-3.0 mg/L, 24 tests

Product #:

LCK390.00

AED Price:

Contact Hach

Hazardous



For the photometric determination of adsorbable organic halogens. Easy to use, packaged in a fully equipped box.

### Specifications

According to standard:	DIN EN ISO 9562
Description:	Adsorbable organic halogens
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.05 - 3.0 mg/L AOX
Method:	Digestion + Iron(III)-Thiocyanate
Number of tests:	24
Parameter:	AOX
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)