



Chloride cuvette test 1-70 mg/L / 70-1000 mg/L Cl, 24 tests

Product #:

LCK311

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	1 - 70 mg/L Cl
Measuring range (2):	70 - 1000 mg/L Cl
Method:	Iron(III)-Thiocyanate
Number of tests:	24
Parameter:	Chloride
Platform:	LCK
Special remarks:	2 measuring ranges with 1 test
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)