



## Sulphide cuvette test 0.1-2.0 mg/L S<sup>2-</sup>, 25 tests

Product #:

LCK653

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	ISO 10530-1991, DIN 38405-D26
Description:	dissolved Sulphide
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.1 - 2.0 mg/L S <sup>2-</sup>
Method:	Dimethyl-p-phenylenediamine
Number of tests:	25
Parameter:	Sulphide
Platform:	LCK
Storage conditions:	15 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)