



Ammonium cuvette test 0.015-2.0 mg/L $\text{NH}_4\text{-N}$, 25 tests

Product #:

LCK304

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017, ISO 23695
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.015 - 2.0 mg/L $\text{NH}_4\text{-N}$
Measuring range (2):	0.02 - 2.5 mg/L NH_4
Method:	Indophenol Blue
Number of tests:	25
Parameter:	Ammonium
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)
Truecal:	Yes

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)