



## Nitrite trace cuvette test 0.0015-0.03 mg/L NO<sub>2</sub>-N, 50 tests

Product #:

LCK541

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

---

### Specifications

According to standard:	EN ISO 26777, DIN 38405 D10
Description:	Nitrite trace cuvette test
Instrument:	DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.0015 - 0.03 mg/L NO <sub>2</sub> -N
Measuring range (2):	0.005 - 0.1 mg/L NO <sub>2</sub>
Method:	Diazotisation
Number of tests:	50
Parameter:	Nitrite, trace
Platform:	LCK
Special remarks:	LCK 541 includes reagent set and LKP341 (50mm semi-micro plastic cuvettes, 10 pcs)
Storage conditions:	2 - 8 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)