



Phosphate (ortho/tot) trace cuvette test 0.01-0.5 mg/L PO₄-P

Product #:

LCS349

AED Price:

Contact Hach

Hazardous



Specifications

According to standard:	ISO 6878-1-1986, DIN 38405 D11-4
Description:	ortho/total Phosphate
Digestion required:	for total P
Measuring range:	0.01 - 0.5 mg/L PO ₄ -P
Measuring range (2):	0.03 - 1.5 mg/L PO ₄
Method:	Phosphormolybdenum Blue
Number of tests:	25
Parameter:	Phosphate, ortho + total
Platform:	LCK
Special remarks:	For the evaluation 50 mm semi-micro cuvettes (LZP341) are needed
Standard method:	EN ISO 6878
Storage conditions:	15 - 25 °C

Required Accessories

- Sample cell, 50 mm rectangular PMMA, 10 pcs. (Item LZP341)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)