



Bühler 1027 Stationary easy wall-mount water sampler

Product #: BU1027.XX.XXXXXX
AED Price: Contact Hach

Basic wall-mount sampler using pressure-vacuum technology. With weatherproof plastic housing, suitable for wall-mounting.

The well organised display and numeric keypad enable the sampler to be rapidly programmed. The instrument can be used for time and volume-proportional sampling. Its simple structure means that it requires very little maintenance. It is weatherproof and can be mounted or fixed on walls.

Pressure-vacuum samplers are ISO 5667 compliant and therefore satisfy the requirements for subsequent reproducible analysis.

The samples can be cooled to 4 °C in a stainless steel refrigerator cabinet (optional).

Ideal for straightforward applications involving the sampling of widely different types of wastewater.

Corrosion-resistant plastic housing

Time and volume proportional composite sampling

Volume accuracy through pressure-vacuum principle

Ideal in combination with stainless steel refrigerator cabinet

Specifications

| | |
|---------------------------|---|
| Cooling system: | External refrigerator option available |
| Housing material: | Reinforced plastic |
| Languages user interface: | HR, SI, GR, GB,NL, IT, DK, SE, PL, ES, RU, DE, FR, PT, CZ, HU, RO, SK, BG, TR |
| Sample container: | 1 x 10 L, 1 x 25 L |
| Sampler type: | Stationary Wall Mount |
| Warranty: | 24 months |