



## Titration Application Pack Water: pH, Alkalinity & Hardness (ISE)

**Product #:** AP0005.AT1222  
**AED Price:** Contact Hach  
**Call for ship date**

Application Pack for the determination of pH, Alkalinity and Calcium & Magnesium Hardness in water. Suited for Titrilab AT1000 Series, Automatic Titrator Model AT1222.

---

### Specifications

Application:	Municipal Applications
Compatible instruments:	AT1222
Parameter:	pH/Alkalinity & Hardness
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.