



## Titration Application Pack Biogas: FOS/TAK

**Product #:** AP0006.AT1102

**AED Price:** Contact Hach

**Call for ship date**

Application Pack for the determination of FOS/TAK (Biogas). Suited for Titrabest AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

---

### Specifications

|                         |   |
|-------------------------|---|
| Application:            | Environmental Applications  |
| Compatible instruments: | AT1102  |
|                         | AT1112  |
|                         | AT1122  |
|                         | AT1222  |
| Parameter:              | FOS/TAK (Biogas)  |
| Type:                   | Application Pack  |
| Warranty:               | 24 months   |
| What's included?:       | Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper. |

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.