



## Titration Application Pack Amperometric: Free & Total Chlorine, Chlorine Dioxide, Chlorite, & Sulfite (Autocat)

**Product #:** AP0007.AT1122  
**AED Price:** Contact Hach  
**Call for ship date**

Application Pack for the determination of free & total Chlorine, Chlorine Dioxide, Chlorite, and Sulphite (Autocat). Suited for Titrilab AT1000 Series, Automatic Titrator Models AT1122/AT1222.

---

### Specifications

Application:	Municipal Applications
Compatible instruments:	AT1122 AT1222
Parameter:	Free & Total Chlorine, Chlorine Dioxide, Chlorite, Sulfite (Autocat)
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

---

### What's included?

Application Pack contains glass syringe, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.