



Formaldehyde trace cuvette test, ISO12460, 0.05 - 3.0 mg/L H_2CO

Product #:

LCS425

AED Price:

Contact Hach

Specifications

According to standard:	ISO12460
Digestion required:	yes
Instrument:	DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.05 - 3.0 mg/L H_2CO
Method:	Acetylacetone
Number of tests:	25
Parameter:	Formaldehyde (trace)
Platform:	LCK
Special remarks:	For the evaluation 50 mm semi-micro rectangle-cuvettes are needed.
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)