



Crucible, Gooch, 40 mL Capacity

| | |
|-------------------|----------------|
| Product #: | 2098641 |
| AED Price: | Contact Hach |
| Available | |

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Specifications

| | |
|-----------|-------|
| Capacity: | 40 mL |
| Quantity: | 1 /pk |