



Beaker, griffin, low form, glass, 2000 mL

| | |
|------------|--------------|
| Product #: | 50054 |
| AED Price: | Contact Hach |
| Available | |

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specifications

| | |
|-------------------|--------------------------|
| Capacity: | 2000 mL |
| Housing material: | Glass |
| Material: | Glass, Griffin, Low Form |
| Units: | each |