



Cylinder, graduated, polymethylpentene, 250 mL \pm 1.4 mL, 2.0 mL divisions

Product #:	217246
AED Price:	Contact Hach
Available	

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated “To Contain/To Deliver” at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Specifications

Capacity:	250 mL
Certifications:	Class B
Housing material:	PMP
Measurement method:	To Deliver/To Contain
Quantity:	1 each
Type:	Single Metric