



## Sample Cell: 1" Square Glass 10 mL matched set of 8

Product #:	2495408
AED Price:	Contact Hach
Available	

Comes as a matched set of 8 cells.

### Use in Combination with Hach Spectrophotometers

Hach spectrophotometers, chemistries, and methods are designed, developed, and manufactured to optimize your results.

---

### Specifications

Footnote:	1 pair of cells is included with purchase of DR 5000 UV-Vis, DR 3800, DR 2800, or DR 2700 Spectrophotometer
Quantity:	8 /pk