



Tubing, C-flex, 6.4 mm (i.d.)

| | |
|-------------------|--------------|
| Product #: | 2327337 |
| AED Price: | Contact Hach |
| Available | |

White, thermoplastic elastomer with excellent chemical resistance and temperature range from -80 to 230°C. Can be sterilized by gamma radiation or autoclave. Supplied in 25-foot sections.

Specifications

| | |
|-------------------|---------|
| Diameter: | 5/16 in |
| Housing material: | C-flex |
| I. D.: | 1/4 in |
| Length: | 25 ft |
| Wall Thickness: | 1/16 in |