



Be Right™



Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Blue Cap/ Stem, 6/pk

Product #:

2928202

AED Price:

Contact Hach

Ships within 3 days

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Specifications

Capacity:	500 mL
Color:	Clear
Material:	LDPE
Quantity:	6 pk
Stem:	Blue Cap/Stem
Unit:	6 pk