



## Test tube/ MPN rack, polyethylene, 20 mm (O.D.), 40 holes

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2497902</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

Forty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

---

### Specifications

|                        |        |
|------------------------|--------|
| Diameter:              | 20 mm  |
| Quantity:              | 1 each |
| Test tube arrangement: | 4 x 10 |
| Test tube holes:       | 40     |