



Pipet, serological, sterile, glass, 5 mL, 25/pk

| | |
|-------------------|----------------|
| Product #: | 2093127 |
| AED Price: | Contact Hach |
| Available | |

Sterile, glass pipet with blue graduation scale. Calibrated “To Deliver”. Complies with ASTM E438, Class B requirements and designed from ASTM Specification E714.

Specifications

| | |
|--------------|--------|
| Capacity: | 5 mL |
| Graduations: | 0.1 mL |
| Quantity: | 25 /pk |