



## Thermometer, Non-Mercury, 0 - 220 °F

<b>Product #:</b>	<b>2676200</b>
AED Price:	Contact Hach
Available	

Compact and durable, non-mercury thermometers enclosed in a protective nylon case. These thermometers are designed for ruggedness and come in many standard ranges. Suitable for all field and portable uses where an unprotected thermometer may break.

---

### Specifications

Quantity:	1 each
Temperature range:	0 °F - 220 °F
Thermometer divisions:	2 °F
Thermometer fluid type:	Non-Mercury