



Test kit, general boiler water, drop count titration, model AL-94 MG-L

Product #:

222403

AED Price:

Contact Hach

Ships within 2 weeks

Hazardous



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus to perform analysis in the field

Uses drop count titration and color disc methods

Specifications

Application:	Boiler Water
Model:	AL-94 MG-L
Parameter:	Drop Count Titration
Platform:	Color Disc