



Test kit, nitrogen , model NI-8

Product #: 224100
AED Price: Contact Hach
Available

Hazardous



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case style:	B
Measuring range:	0 - 3.0 mg/L NH ₃ -N
Method:	Color disc/Nessler reagent
Number of tests:	100
Parameter:	Nitrogen, Ammonia, mid range - For freshwater
Platform:	Color Disc
Smallest increment:	0.1
Test kit model:	NI-8