



## Filter, acrocap, pore size 0.2 µm

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2323810</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

---

### Specifications

|              |                 |
|--------------|-----------------|
| Diameter:    | NA              |
| Filter Type: | Filter, acrocap |
| Packing:     | NA              |
| Pore size:   | 0.2 µm          |