



## Bottle rack, plastic-coated, 20 compartments

<b>Product #:</b>	<b>2551801</b>
AED Price:	Contact Hach
Available	

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

---

### Specifications

Quantity:	1
Special Features:	Rack, 20 Bottle Compartments