



pH Buffer Solution, pH 4.01, 500 mL

Product #:	2283449.99
AED Price:	Contact Hach
Available	

Quality Buffer Solution, colour-coded red for easy identification. NIST traceable. pH is 4.01 ± 0.02 at 25°C.

Specifications

Color coded?:	Yes
Concentration:	4.01 ± 0.02 pH
Description:	Quality Buffer Solution, colour-coded red (NIST)
Package Type:	Bottle
Parameter:	pH
pH value(s):	4.01
Platform:	Solution
Quantity:	500 mL