



## Dissolved Iron Chemkey Reagents (box of 25)

**Product #:** 8636000  
**AED Price:** Contact Hach

**Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.**

Dissolved Iron Chemkeys, range 0.05 - 3.00 mg/L Fe. For SL1000 Portable Parallel Analyser (PPA).

**Current SL1000 instrument customers who want to test Dissolved Iron need to perform a software update that can be located on the instrument download section.**

---

### Specifications

Measuring range:	0.05 - 3.00 mg/L Fe
Parameter:	Dissolved Iron
Platform:	Chemkeys
Quantity:	pk/25
What's included?:	25 Chemkeys

---

### What's included?

25 Chemkeys