



Low Range Total Alkalinity Chemkey Reagents (box of 25)

Product #: 8636200
AED Price: Contact Hach

Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.

Low Range Total Alkalinity Chemkeys, range 20 - 200 mg/L CaCO_3 . For SL1000 Portable Parallel Analyser (PPA).

Current SL1000 instrument customers who want to test Alkalinity need to perform a software update that can be located on the instrument download section.

Specifications

Measuring range:	20 - 200 mg/L CaCO_3
Method:	Colorimetric
Method number:	10280
Parameter:	Alkalinity, LR
Platform:	Chemkeys
Quantity:	pk/25