



Cuvette Test Permanganate Index 0.5 - 10 mg/L O₂ (CODMn), 25 tests

Product #:

LCK394

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Specifications

According to standard:	ISO 8467
Description:	Permanganate Index (CODMn)
Digestion required:	Yes
EPA compliant:	N/A
Instrument:	DR1900, DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.5 - 10 mg/L O ₂
Method:	Permanganate
Number of tests:	25
Parameter:	CODMn
Platform:	LCK
Standard method:	ISO 8467
Storage conditions:	15 - 25 °C (protect from light)