



Manganese Reagent Set 0.005 - 0.7 mg/L Mn

Product #:

LCW632

Hazardous

AED Price:

Contact Hach



Cyanide-free reagent set for determination of Manganese with TMB

Specifications

Instrument:	DR3900, DR6000, DR2800, DR3800, DR5000
Measuring range:	0.005 - 0.7 mg/L Mn
Measuring range (2):	Ultra Low Range
Method:	TMB
Number of tests:	50
Parameter:	Manganese
Platform:	LCW
Special remarks:	50 mm rectangle-cuvette LZP341 needed
Storage conditions:	4 - 8 °C (protect from light)