



Peracetic Acid Chemkey Reagents (box of 25)

Product #: 8635200
AED Price: Contact Hach

Peracetic Acid Chemkeys, range 0.04 - 50.0 mg/L PAA. For SL1000 Portable Parallel Analyzer (PPA).

Specifications

Instrument:	SL1000
Measuring range:	0.04 - 50.0 mg/L PAA
Method number:	10297
Parameter:	Peracetic Acid
Platform:	Chemkeys
Quantity:	pk/25
What's included?:	25 Chemkeys

What's included?

25 Chemkeys